

## Product Datasheet

# Monkeypox virus/MPXV E8L Recombinant Protein (orb1671357)

**Description** Monkeypox virus/MPXV E8L Recombinant Protein

**Tested Applications** ELISA, SDS-PAGE, WB

**Preservatives** Lyophilized from PBS pH 7.4, 1mM EDTA, 4%trehalose, 1% mannitol. Reconstitute in sterile water for a stock solution.

**Form/Appearance** Lyophilized

**Storage** Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

**Tag** C-His

**Note** For research use only

**Protein Sequence** Pro2-Lys274

**Purity** > 90% as determined by SDS-PAGE


**Source** Mammalian Cells

**MW** 34.94 kDa


**Uniprot ID** [Q8V4Y0](#)

**NCBI** [Q8V4Y0](#)

**Expiration Date** 6 months from date of receipt.



Anti-MPXV E8L Antibody binds with E8L pr...



SDS-PAGE