

Product Datasheet

Monkeypox virus/MPXV B5R Recombinant Protein (orb1671361)

Description	Monkeypox virus/MPXV B5R Recombinant Protein
Tested Applications	ELISA, SDS-PAGE, WB
Preservatives	Lyophilized from PBS pH 7.4 [0.02% NLS, 1mM EDTA, 4% trehalose [1% mannitol. Reconstitute in sterile water for a stock solution.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Tag	N-His
Note	For research use only
Protein Sequence	Asn26-Thr371
Purity	> 90% as determined by SDS-PAGE
Source	E. coli
MW	42.71 kDa
Uniprot ID	Q8V4S3
NCBI	Q8V4S3
Expiration Date	6 months from date of receipt.