



Product Datasheet

Monkeypox virus/MPXV A29L Recombinant Protein (orb1671365)

Description	Monkeypox virus/MPXV A29L Recombinant Protein
Tested Applications	ELISA, SDS-PAGE, WB
Preservatives	Lyophilized from PBS pH 7.4, 1mM EDTA, 4%trehalose, 1% mannitol. Reconstitute in sterile water for a stock solution.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Tag	C-His
Note	For research use only
Protein Sequence	Asp2-Glu110
Purity	> 90% as determined by SDS-PAGE
Source	Mammalian Cells
MW	15.33 kDa
Uniprot ID	Q77HM6
NCBI	Q77HM6
Expiration Date	6 months from date of receipt.



Anti-MPXV A29L Antibody binds with A29L ...



SDS-PAGE