

Product Datasheet

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody [T3P1D2] (orb1671389)

Description	SARS-CoV/SARS-CoV-2 Variant (Delta,	
Species/Host	Antibody [T3P1D2]	
Reactivity	Camelus	
Conjugation	Virus	
Tested Applications	Unconjugated	
Immunogen	SARS-CoV-2 variant Spike proteins.	ELISA Validation with Trimers of SARS1 a...
Target	S	
Preservatives	SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody is supplied in PBS.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Isotype	sdAb	
Clonality	Recombinant	
Uniprot ID	P0DTC2	
NCBI	QHD43419	
Expiration Date	12 months from date of receipt.	