



Product Datasheet

GAG Antibody (Biotin) (orb1671393)

Description	GAG Antibody (Biotin)	
Species/Host	Mouse	
Reactivity	Virus	
Conjugation	Biotin	
Tested Applications	ELISA, WB	
Immunogen	Mouse monoclonal HIV-1 p24 antibody was raised against a recombinant full-length HIV-1 p24 protein.	
Target	GAG	Western Blot Validation with HIV-1 p24 P... 
Preservatives	HIV-1 p24 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	HIV-1 p24 Monoclonal Antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	ELISA Validation of HIV-1 p24 Antibody w...
Application notes	WB: 1 µg/mL. Anti-HIV-1 p24 antibody detects HIV-1 p24 protein by ELISA and WB. All other applications and species not yet tested.	
Isotype	IgG1	
Clonality	Monoclonal	
MW	24 kDa	
Uniprot ID	P04591	
NCBI	AAB50258	
Expiration Date	12 months from date of receipt.	