


---


## Product Datasheet

### GAG Antibody (orb1671394)

<b>Description</b>	GAG Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Mouse monoclonal HIV-1 p24 antibody was raised against a recombinant full-length HIV-1 p24 protein.
<b>Target</b>	GAG
<b>Preservatives</b>	HIV-1 p24 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	HIV-1 p24 Monoclonal Antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: 1 µg/mL. Anti-HIV-1 p24 antibody detects HIV-1 p24 protein by ELISA and WB. All other applications and species not yet tested.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	24 kDa
<b>Uniprot ID</b>	<a href="#">P04591</a>
<b>NCBI</b>	<a href="#">AAB50258</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot Validation with HIV-1 p24 P...



ELISA Validation of HIV-1 p24 Antibody w...