


---

## Product Datasheet

### GAG Antibody (HRP) (orb1671395)

Description	GAG Antibody (HRP)	KDa
<b>Species/Host</b>	Mouse	250 -
<b>Reactivity</b>	Virus	72 -
<b>Conjugation</b>	HRP	55 -
<b>Tested Applications</b>	ELISA, WB	36 -
<b>Immunogen</b>	Mouse monoclonal HIV-1 p24 antibody was raised against a recombinant full-length HIV-1 p24 protein.	28 -
<b>Target</b>	GAG	17 -
<b>Preservatives</b>	HIV-1 p24 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.	←
<b>Form/Appearance</b>	Liquid	Western Blot Validation with HIV-1 p24 P...
<b>Concentration</b>	1 mg/mL	
<b>Storage</b>	HIV-1 p24 Monoclonal Antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	ELISA Validation of HIV-1 p24 Antibody w...
<b>Note</b>	For research use only	
<b>Application notes</b>	WB: 1 µg/mL. Anti-HIV-1 p24 antibody detects HIV-1 p24 protein by ELISA and WB. All other applications and species not yet tested.	
<b>Isotype</b>	IgG1	
<b>Clonality</b>	Monoclonal	
<b>MW</b>	24 kDa	
<b>Uniprot ID</b>	<a href="#">P04591</a>	
<b>NCBI</b>	<a href="#">AAB50258</a>	
<b>Expiration Date</b>	12 months from date of receipt.	