


Product Datasheet

Il2ra Antibody (orb1671485)


Description	Il2ra Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Depletion, FC
Immunogen	The original version of this antibody was raised to recognize mouse CD25.
Target	Il2ra
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P01590
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-CD25 antibody...



Western Blot using Anti-CD25 antibody 2E...



Immunofluorescence staining of mouse spl...