

Product Datasheet

Spike protein (RBD) Antibody [Sb#14], Human IgG1-Fc Fusion (orb1671513)

Description	Spike protein (RBD) Antibody [Sb#14], Human IgG1-Fc Fusion
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FA
Immunogen	This clone was originally isolated in a form of a synthetic nanobody (sybody) via a 'target swap' selection procedure against RBD-vYFP using ribosomal display and against RBD-Fc fusion during phage display rounds.
Target	S
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Tag	Fc Silent
Note	For research use only
Isotype	IgG1-Fc Fusion
Clonality	Recombinant
Uniprot ID	PODTC2
Expiration Date	12 months from date of receipt.