

## Product Datasheet

# Spike protein (RBD) Antibody [Sb#14], Mouse IgG1-Fc Fusion (orb1671514)

<b>Description</b>	Spike protein (RBD) Antibody [Sb#14], Mouse IgG1-Fc Fusion
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA
<b>Immunogen</b>	This clone was originally isolated in a form of a synthetic nanobody (sybody) via a 'target swap' selection procedure against RBD-vYFP using ribosomal display and against RBD-Fc fusion during phage display rounds.
<b>Target</b>	S
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1-Fc Fusion
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>PODTC2</b>
<b>Expiration Date</b>	12 months from date of receipt.