


Product Datasheet

CD38 Antibody (orb1671515)

Description	CD38 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	The original antibody was generated by immunizing CAF1 mice with human thymocytes.
Target	CD38
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P28907
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed Dau...