



Product Datasheet

FOXP2 Antibody (orb1671517)

Description	FOXP2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing FOXP2 protein under non-denaturing conditions.
Target	FOXP2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	O15409
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeL...



Flow cytometry using the anti-FOXP2 anti...