



Product Datasheet

FOXP4 Antibody (orb1671521)

Description	FOXP4 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing FOXP4 protein under non-denaturing conditions.
Target	FOXP4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q81VH2
Expiration Date	12 months from date of receipt.



Flow cytometry using the anti-FOXP4 anti...



Immunofluorescence staining of fixed HeL...



Western Blot using anti-FOXP4 antibody N...