



Product Datasheet

SFMBT1 Antibody (orb1671527)

Description	SFMBT1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing SFMBT1 protein under non-denaturing conditions.
Target	SFMBT1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9UHJ3
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of A431 cell...



Flow cytometry using the Anti-SFMBT1 ant...