



Product Datasheet

PRDM9 Antibody (orb1671533)


Description	PRDM9 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing PRDM9 protein under non-denaturing conditions.
Target	PRDM9
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9NQV7
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Immunofluorescence staining of HEK293 ce...



Flow cytometry using the Anti-PRDM9 anti...