



Product Datasheet

CBX1 Antibody (orb1671545)


Description	CBX1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing CBX1 protein under non-denaturing conditions.
Target	CBX1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P83916
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeL...



Western Blot using anti-CBX1 antibody R ...



Flow cytometry using the anti-CBX1 antib...