


Product Datasheet

Covid-19 & SARS-CoV Nucleoprotein Antibody [CR3018 (03-018)], Rabbit IgG (orb1671555)

Description	Covid-19 & SARS-CoV Nucleoprotein Antibody [CR3018 (03-018)], Rabbit IgG
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	The original antibody was generated by cloning the variable regions of the scFvs selected from phage display libraries into separate vectors for IgG1 heavy-chain and light-chain expression. The harvested supernatants were then purified on protein A columns. The original immunogen was the whole irradiated virion.
Target	N
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P0DTC9, P59595
Expiration Date	12 months from date of receipt.



Binding curve of anti-SARS-Cov-2 (COVID-...