



## Product Datasheet

### TNFRSF9 Antibody (orb1671559)

<b>Description</b>	TNFRSF9 Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, ELISA, FC, IHC
<b>Immunogen</b>	The original antibody was generated by immunizing BALB/c mice with a fusion protein of human 4-1BB antigen and glutathione s transferase (GST).
<b>Target</b>	TNFRSF9
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q07011</a>
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence staining of fixed Jur...



Western Blot using anti-4-1BB antibody 4...