



Product Datasheet

CD80 Antibody (orb1671561)

Description	CD80 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA
Immunogen	This antibody was generated by immunization of cynomolgus monkeys with human B7.1 antigen.
Target	CD80
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	P33681
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed Dau...



Flow cytometry using the anti-CD80 antib...