



## Product Datasheet

### CD74 Antibody (orb1671568)

<b>Description</b>	CD74 Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IP, RIA, WB
<b>Immunogen</b>	This antibody was generated by immunizing female BALB/c mice with whole Raji cells or a crude membrane preparation.
<b>Target</b>	CD74
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P04233</b>
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot using anti-CD74 antibody LL...



Immunofluorescence staining of Daudi cel...