

Product Datasheet

COVID-19 & SARS-CoV S glycoprotein Antibody [CR3022], Mouse IgG2b (orb1671574)

Description	COVID-19 & SARS-CoV S glycoprotein Antibody [CR3022], Mouse IgG2b
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	Crystallography, ELISA, IF
Immunogen	The original monoclonal antibody was generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV.
Target	S
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b kappa
Clonality	Recombinant
Uniprot ID	P59594
Expiration Date	12 months from date of receipt.