


Product Datasheet

CD40LG Biosimilar Antibody (orb1671579)

Description	CD40LG Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, IF
Immunogen	The original monoclonal antibody was generated by immunization with a soluble fusion protein of human gp39 (gp39-CD8).
Target	CD40LG
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P29965
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed Jurkat cells.