



Product Datasheet

CXCR2 Antibody (orb1671589)

Description	CXCR2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	This antibody was raised by immunizing mice with human IL-8 (interleukin-8) receptor, type B.
Target	CXCR2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P25025
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-IL-8RB ant...



Immunofluorescence staining of fixed hum...