



Product Datasheet

Ccr9 Antibody (orb1671594)

Description	Ccr9 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, NeA
Immunogen	This antibody was raised by immunizing rats with KLH-conjugated linear peptide corresponding to the N-terminus of mouse CCR9.
Target	Ccr9
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9WUT7
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-CCR9 antibody...



Immunofluorescence staining of fixed RAW...