



Product Datasheet

CCR5 Antibody (orb1671595)

Description	CCR5 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	This antibody was raised by immunizing mice with CCR5.
Target	CCR5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P51681
Expiration Date	12 months from date of receipt.




Immunofluorescence staining of fixed U93...



Western Blot using anti-CCR5 (PKC site)

...



Flow cytometry using the Anti-CCR5 (PKC ...)