



Product Datasheet

CCR5 Antibody (orb1671601)

Description	CCR5 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IP, WB
Immunogen	This antibody was raised by immunizing mice with intact RBL-CCR5 cells.
Target	CCR5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P51681
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-CCR5 antibody...



Immunofluorescence staining of fixed U93...