



Product Datasheet

C5AR1 Antibody (orb1671603)

Description	C5AR1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	This antibody was raised by immunizing mice a synthetic C5aR C-terminal phosphopeptide that contained phosphoserine instead of serine at position 334 and also at the neighboring phosphorylation sites Ser-332 and Ser-338.
Target	C5AR1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P21730
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed U937 cells.



Flow cytometry using the Anti-C5aR antibody.