



Product Datasheet

VCAM1 Antibody (orb1671611)

Description	VCAM1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunizing Balb/c mice against activated human umbilical vein endothelial cells (HUVEC).
Target	VCAM1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P19320
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed hum...



Flow cytometry using the Anti-CD106 anti...