



Product Datasheet

CD81 Antibody (orb1671615)

Description	CD81 Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	This antibody was raised by immunizing mice with a B-Cell line derived from a Burkitt lymphoma.
Target	CD81
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P60033
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-CD81 antibody.



Immunofluorescence staining of fixed cells.