


Product Datasheet

CD63 Antibody (orb1671621)

| | |
|----------------------------|---|
| Description | CD63 Antibody |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | FC, IF, IHC, IP, WB |
| Immunogen | This antibody was raised by immunizing BALB/c mice with a plasma membrane preparation of the human melanoma MeWo cells. |
| Target | CD63 |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | batch dependent |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Isotype | IgG kappa |
| Clonality | Recombinant |
| Uniprot ID | P08962 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence staining of fixed A43...



Western Blot using anti-CD63 antibody NK...