



Product Datasheet

CD37 Antibody (orb1671630)

Description	CD37 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Immunogen	This antibody was raised by immunizing mice with spleen cells of a patient with hairy cell leukemia.
Target	CD37
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P11049
Expiration Date	12 months from date of receipt.



Western Blot using anti-CD37 antibody IP...



Immunofluorescence staining of fixed Dau...