



Product Datasheet

CD1A Antibody (orb1671633)

Description	CD1A Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP
Immunogen	This antibody was raised by immunizing mice with the human CD1a protein.
Target	CD1A
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P06126
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed Jurkat cells.



Flow cytometry using the anti-CD1a antibody.