



Product Datasheet

CD19 Antibody (orb1671645)

Description	CD19 Antibody
Reactivity	Human
Conjugation	Unconjugated
Immunogen	This antibody was raised by immunising mice against CD19.
Target	CD19
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P15391
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-CD19 antibody.



Immunofluorescence staining of fixed human cells.