



Product Datasheet

PTPRC Antibody (orb1671648)

Description	PTPRC Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	This antibody was raised by immunising mice with an EBV-transformed B cell line from a patient with MHC class II deficiency.
Target	PTPRC
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P08575
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of Jukat cells...



Flow cytometry using the Anti-CD45RA ant...