



Product Datasheet

TSPAN1 Antibody (orb1671657)

Description	TSPAN1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	This antibody was raised by immunising BALB/c mice with the major extracellular domain of tetraspanin 1.
Target	TSPAN1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	O60635
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed A43...



Flow cytometry using the Anti-Tetraspani...