


Product Datasheet

NCAM1 Antibody (orb1671664)

Description	NCAM1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	This antibody was raised against human NCAM/CD56
Target	NCAM1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P13591
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of Kelly cel...



Western Blot using anti-CD56 antibody NC...