

Product Datasheet

Lewis a,b Antibody [HEA164], Rabbit IgG (orb1671665)

Description Lewis a,b Antibody [HEA164], Rabbit IgG

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, IHC, IP, WB

Immunogen This antibody was raised by immunising mice with human carcinoma cells.

Preservatives PBS with 0.02% Proclin 300.

Concentration batch dependent

Storage Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.

Note For research use only


Isotype IgG kappa

Clonality Recombinant

Expiration Date 12 months from date of receipt.



Western Blot using anti-Lewis a-b antibody



Immunofluorescence staining of fixed Caco-2 cells