



Product Datasheet

FCER1G Antibody (orb1671675)

Description	FCER1G Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant human protein, N-terminus.
Target	FCER1G
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P30273
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of U937 cell...



Flow-cytometry using the Anti-Fc epsilon...