


Product Datasheet

Cd226 Antibody (orb1671694)

Description	Cd226 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	This antibody was raised by immunising Lewis and Lou/M rats with the murine Th1 cell clone AE7 and Th1 cell lines derived from 5B6 and DO11.10 TCR transgenic mice, before screening for reactivity with Th1 clones and a lack of reactivity to Th2/Th0 cell lines.
Target	Cd226
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q8K4F0
Expiration Date	12 months from date of receipt.



Western Blot using anti-CD266 antibody 1...



Flow cytometry using the Anti-CD226 anti...