


## Product Datasheet

### Pvr Antibody (orb1671696)


<b>Description</b>	Pvr Antibody
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC
<b>Immunogen</b>	This antibody was raised by immunising rats with recombinant murine CD155, consisting of the entire ectodomain expressed as a HIS-tagged protein.
<b>Target</b>	Pvr
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q8K094</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot using anti-CD155 antibody 3...



Immunofluorescence staining of RAW 264.7...



Binding curve of anti-CD155 antibody 3F1...