



Product Datasheet

Lag3 Antibody (orb1671700)

Description	Lag3 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IP
Immunogen	This antibody was raised by immunising Lewis rats with a murine LAG-3:Ig fusion protein.
Target	Lag3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q61790
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of mouse spl...



Binding curve of anti-Lag3 antibody C9B7...