



Product Datasheet

Slc3a2 Antibody (orb1671702)

Description	Slc3a2 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP, Stim
Immunogen	This antibody was raised by immunising rats with an IL 2-producing subline of the mouse T lymphoma EL4.
Target	Slc3a2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P10852
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of mouse spl...



Western Blot using anti-CD98 antibody RL...