



## Product Datasheet

### IL9R Antibody (orb1671707)

<b>Description</b>	IL9R Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, IP
<b>Immunogen</b>	This antibody was raised by immunising DBA/2 mice with P815 mastocytoma cells transfected with human IL-9R cDNA.
<b>Target</b>	IL9R
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q01113</a>
<b>Expiration Date</b>	12 months from date of receipt.



Flow-cytometry using the Anti-IL-9R anti...



Immunofluorescence staining of Jurkat ce...