



Product Datasheet

ITGAL Antibody (orb1671720)

Description	ITGAL Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC
Immunogen	Spleen cells from BALB/c mice immunised with EBV-transformed human lymphoblastoid cells were fused with P3-NSUIAg4.1 myeloma cells and the resulting hybridomas were screened for antibody production by a radioimmunoassay on T and B cell lines.
Target	ITGAL
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P20701
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of Jurkat ce...



Flow cytometry using the Anti-CD11a anti...