



Product Datasheet

Tfrc Antibody (orb1671722)

Description	Tfrc Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Depletion, FC
Immunogen	The antibody was obtained from a fusion of the murine myeloma S 194/5.XXO.BU. 1 with spleen cells from a BDIX rat immunized with the Friend erythroleukemia 745.6.
Target	Tfrc
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q62351
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of RAW 264.7...



Flow-cytometry using the Anti-Transferrin...