



Product Datasheet

RHD Antibody (orb1671729)

Description	RHD Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Cytotoxic Assay, FC
Immunogen	A Rh- volunteer blood donor (r/r) was injected with 5ml of Rh+ erythrocytes (R1/r) and the peripheral lymphocytes were isolated and transformed with 1 ml of culture supernatant of the EBV marmoset B96-8 cell line.
Target	RHD
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	Q02161
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed K562...



Flow cytometry using the Anti-Rhesus D a...