



Product Datasheet

ITGB2 Antibody (orb1671731)

Description	ITGB2 Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IP, RIA
Immunogen	Blood was obtained from 4 healthy donors and the mononuclear cells harvested and grown. A BALB/c mouse was injected with 10^6 monocytes three times at monthly intervals, and four times with 5×10^6 monocytes on each of the four days before fusion. The spleen cells were then fused to P3U1 myeloma cells to generate hybridomas and RIA was used to identify lines of interest.
Target	ITGB2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P05107
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-CD18 antibody...



Immunofluorescence staining of human per...