


Product Datasheet

TAG-72 Antibody [Minretumomab (CC49)], Rabbit IgG - Research Grade Biosimilar (orb1671741)

Description	TAG-72 Antibody [Minretumomab (CC49)], Rabbit IgG - Research Grade Biosimilar
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, RIA, WB
Immunogen	This antibody was raised by immunising BALB/c mice with TAG-72 purified from LS-174T human carcinoma xenograft.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of Jurkat ce...