



## Product Datasheet

### Ms4a1 Antibody (orb1671745)

<b>Description</b>	Ms4a1 Antibody
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Depletion, FA, FC, IF, WB
<b>Immunogen</b>	This antibody was generated by immunizing CD20 <sup>-/-</sup> mice with CD20-transfected cell lines (300.18 and 70Z3) and CD20-peptide-keyhole limpet hemocyanin conjugate.
<b>Target</b>	Ms4a1
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P19437</b>
<b>Expiration Date</b>	12 months from date of receipt.



Flow-cytometry using the Anti-CD20 antib...



Immunofluorescence staining of U937 cell...