



Product Datasheet

Tnfrsf18 Antibody (orb1671762)

Description	Tnfrsf18 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, Functional Studies, IP, WB
Immunogen	This antibody was raised by immunising Wistar rats with the CD25+CD4+ T cell line. Hybridoma cells were generated by fusing spleen cells from the immunized rat with P3X63Ag8.653 myeloma cells, assessed for secreting antibodies capable of neutralizing in vitro CD25+CD4+ T cell-mediated suppression of anti-CD3-stimulated T cell proliferation ⁹ , and then cloned. This DTA-1 clone secreted a monoclonal rat immunoglobulin G2a (IgG2a) that showed high neutralizing activity.
Target	Tnfrsf18
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	O35714
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed mou...



Flow cytometry using the Anti-GITR antib...